



Biomed Healthtech Pvt. Ltd.

Mumbai, Ahmedabad

BINOCULAR LOUPES



Optional clip-in correction frame



Professional Headband

HR 2.5x Specifications			
Working distance in mm	340	420	520
Field of view in mm	90	130	180
Depth of field in mm	110	180	260

HRP Specifications			
Working distance in mm	HRP 3.5x/420	HRP 4x/340	HRP 6x/340
Field of view in mm	65	50	40
Depth of field in mm	60	40	30

- Various Magnifications **2.5X, 3.5X, 4X, 6X**
- Multi-Coated Achromatic optics for crisp, bright images with colour-correction
- High magnification with a large field of view and excellent depth of field.
- Independent adjustment of left and right optics for perfect PD adjustment
- Lightweight
- Water-resistant.
- **Heine Binocular Loupe can be attached with Loupe Light2**



- Compact and Light Weight LED Loupe Light
- Coaxial Illumination
- Offers Completely Homogenous Illumination
- Illumination – 55000 Lux
- LED Operating Life – 50000 hrs
- Operating Time – 9 hr with rechargeable battery pack

[VIDEO](#)

[DETAILS](#)

[EMAIL](#)

[WhatsApp](#)

[WEBSITE](#)