



DERMATOSCOPES

DELTA 30

DELTA one



[VIDEO](#)

[VIDEO](#)

[DETAILS](#)

[EMAIL](#)

[WhatsApp](#)

[WEBSITE](#)



Biomed Healthtech Pvt. Ltd.
Mumbai, Ahmedabad

FEATURES	DELTA 30	DELTA one
Polarisation	Polarised & Non-Polarised With Single Switch	Polarised & Non-Polarised With Single Switch
Illumination	LED HQ	LED HQ
Free Field Of View (Diameter In Contact Mode)	30 Mm	21 Mm
Magnification	10 X	10 X
Diaptor Range	-4 to +4 dpt	-4 to +4 dpt
Brightness Control	3 Level of Light Intensity	3 Level of Light Intensity
Working Hours	Min Brightness : 640 min. Max Brightness : 160 min.	Min Brightness : 460 min. Max Brightness : 80 min.
Full Charging Time	120 min.	75 min.
Handle Options / Power Sources	<ul style="list-style-type: none"> • USB Cable • Table Charger 	<ul style="list-style-type: none"> • USB Cable
Contact Plate	Contact Plate With Scale	Contact Plate With Scale
Digital Documentation	With adaptor for Smart Phone and HEINE DERM app	With adaptor for Smart Phone and HEINE DERM app

Digital Documentation

HEINE DERM app

- Storage and Management of Dermoscopic Digital Images.
- Simple Assignment of the Images Using a Bodymap.
- Annotation with Comments, Risk Levels, Diagnoses and Treatments.
- Intuitive Workflow Due To the Clear User Interface.
- Easily Compare Images From Past Examinations.
- Sharing function

Optional Accessories:

Smart Phone Connector

Android

iPhone



Small Contact Plate



[DETAILS](#)

[EMAIL](#)

[WhatsApp](#)

[WEBSITE](#)