




**Biomed Healthtech Pvt. Ltd.**

**Mumbai, Ahmedabad**

## ML4 LED HEADLIGHT



Comparison of homogeneity  
[01] Competitor product  
[02] HEINE ML4 LED HeadLight

- LED In HEINE Quality
- Comfortable And Secure Fit
- Shadow Free
- Illumination Typ. 90,000 Lux
- Adjustable Illumination Spot Size From 30 mm To 80 mm Spot Size Range at 420 mm Working Distance
- Bright And Homogenous Illumination
- Stepless Light Intensity Control
- Coaxial Illumination
- The Professional L Headband Has Multiple Adjustment Points And Soft Padding For Comfort And Stability
- Optional Polarisation Filter
- Heine Binocular Loupes (Magnifications: 2.5X, 3.5X, 4X, 6X) can be attached.

### Power Source:

- Headband Mounted Battery (Cordless)(Operating time 3.5 hours)
- mpack Rechargeable Battery Transformer (Operating time Typ.8.5 hours)
- Plug-In Transformer

[VIDEO](#)

[DETAILS](#)

[EMAIL](#)

[WhatsApp](#)

[WEBSITE](#)