




**Biomed Healthtech Pvt. Ltd.**

**Mumbai, Ahmedabad**

**DIRECT OPHTHALMOSCOPE- BETA 200**

With **XHL & LED HQ** Illumination



- Aspherical Optical System
- Aluminium Frame System
- Eliminates Corneal Reflex
- Ideal For Small Pupils
- Range Of Lenses:
  - 27 Lenses From -35D To 40D
  - + In 1 D Steps: 1 – 10 15 20
  - In 1 D Steps: 1 – 10 15 20 25 35D

**Power Sources Options:**



**Battery Handle**



**USB Rechargeable Handle & Charger**



**Wall Mounted Unit  
(Electrically Operated)**

**STREAK RETINOSCOPE - BETA 200**

With **XHL & LED HQ** Illumination



**Optional Accessories:  
Dynamic Retinoscopy  
Cards**



- Parastop
- Highly Visible, Bright Fundus Reflex
- Metal Controls
- Single Control For Vergence & Rotation
- Integrated Polarisation Filter
- Detachable Brow Rest

**COMBINED DIAGNOSTIC SET**



[DETAILS](#)

[EMAIL](#)

[WhatsApp](#)

[WEBSITE](#)