



Biomed Healthtech Pvt. Ltd.

Mumbai, Ahmedabad

OTOSCOPE F.O. - BETA 400 LED



- Exclusive **4.2 X Magnification** Offers more details of Anatomy
- Fiber Optic Illumination Ensures Homogeneous, Very Bright Illumination and an Unobstructed View Of The Ear Canal and Tympanum
- All-Metal Housing
- Integrated Insufflation Port
- Large Field Of View Gives The Necessary Overview And Perfect Orientation During Examinations
- Multi-Coated Precision Optics Offer High Resolution and Distortion Free Image
- No Need To Replace Bulbs

Optional Accessories:

**Disposable Tips 2.5 mm / 4 mm
in a Pack of 250 / 1000 Tips**



Pneumatic Bulb



[VIDEO](#)

[DETAILS](#)

[EMAIL](#)

[WhatsApp](#)

[WEBSITE](#)